## **ULTRASOUND THERAPY**

## HS-502

HS-502 uses 1MHz of ultrasonic waves as treatment frequency for relieving pain. It applies constant supersonic waves to the treatment part, providing the maximum effect for treatment with crystal vibrator. Also, it is very easy to operate and control the equipment and it provides various pulse modes.

## **Specifications**

Power	AC 220V, 50/60Hz
Power Consumption	75VA
Safety, Isolation	Class I, Type BF
Output Frequency	1Mhz
Output Power	0.3~1.5W/cm²
Effective Radiated Area	3.5cm²
Timer	1~30min
Output Mode	5(5 sec on/off), 10(2 sec on/off), 30(5 sec on/2 sec off), continuous
Dimension	250x280x110(H)mm
Weight	5Kg

## Accessory



Hand piece

