

ORTHOPEDICS TRACTION UNIT

HT-101 Biotrack, HT-102 Biotrack

The electrically-powered orthopedics traction unit is designed to use softly and easily. This appliance can track the cervical and the waist vertebrae individually at the same time. HT-101 is for one patient and HT-102 for two patients. An Emergency stop switch for a patient is provided.

Feature

- 4 traction modes · Computerized control.
- Hot-pack effect in the back area of the treatment table.
- Convenient operation with shift lock.
- If a patient moves during treatment, the device adjusts automatically to the Patient's position.
- It is always maintaining fixed pulling power for treatment.


Specifications

Power	AC 220V, 60Hz
Power Consumption	120W(HT-101)/200W(HT-102)
Safety, Isolation	Class I, Type B
Pulling Power	1~80Kg
Pulling Speed	7mm/sec
Timer	1~99min
Hold/Standby Time	1~99sec
Thermal effect	[Level 1] 40°C / [Level 2] 60°C
Dimension	2425x580x2640(H)mm
Weight	60kg(Control Unit), 65kg(Bed)



* HT-102S (with one bed), HT-102D (with two beds)



HT-101
Single system



HT-102
Dual system



HT-101 Biotrack



HT-101 Biotrack
Single Control Panel



HT-102 Biotrack
Dual Control Panel



HT-102 Biotrack

Accessory

